



 Operating Instructions

 Mode d'emploi

 Bedienungsanleitung

 Bruksanvisning

 Handleiding voor de Bediening

 Bruksanvisning

 Bruksanvisning

 Instrucciones de funcionamiento

 Instruções para Utilização

 Istruzioni per l'uso

 Οδηγίες Λειτουργίας



*2100
Classic*

Clinical Autoclaves



P3 - Decals, displays and controls
P4 - Introduction
- Key to Pictures
P5 - Getting Started
P6 - Continued Operation
P7 - Troubleshooting
P8 - Technical Specification
P9 - Additional Information



S3 - Dekalkeringer, displayer og kontroller
S40 - Introduksjon
- Bildeforklaring
S41 - Igangsetting
S42 - Fortsatt bruk
S43 - Feilsøking
S44 - Teknisk spesifisering
S45 - Ytterligere informasjon



P3 - Autocollants, panneaux et commandes
P10 - Introduction
- Légendes des images
P11 - Prise en main
P12 - Fonctionnement continu
P13 - Dépistage des anomalies
P14 - Caractéristiques techniques
P15 - Renseignements complémentaires



P3 - Calcomanías, visualizadores y controles
P46 - Introducción
- Clave a las fotografías
P47 - Preparación para el funcionamiento
P48 - Funcionamiento continuo
P49 - Localización de averías
P50 - Especificaciones técnicas
P51 - Información adicional



S.3 - Symbole, Anzeigen und Bedienungselemente
S.16 - Einleitung
- Zeichenerklärung zu Bildern
S.17 - Inbetriebnahme
S.18 - Kontinuierlicher Betrieb
S.19 - Fehlererkennung und beseitigung
S.20 - Technische Daten
S.21 - Weitere Informationen



P3 - Dísticos, mostradores e comandos
P52 - Introdução
- Legendas das imagens
P53 - Como começar
P54 - Funcionamento prolongado
P55 - Resolução de problemas
P56 - Especificações Técnicas
P57 - Informação Adicional



S3 - Mærkater, displays og betjeningsanordninger
S22 - Indledning
- Nøgle til billeder
S23 - Sådan kommer man i gang
S24 - Fortsat drift
S25 - Fejlfinding
S26 - Tekniske specifikationer
S27 - Supplerende informationer



P3 - Decalcomanie, display e comandi
P58 - Introduzione
- Legenda delle Illustrazioni
P59 - Procedura di Avviamento
P60 - Funzionamento Ininterrotto
P61 - Ricerca dei Guasti
P62 - Specifica Tecnica
P63 - Ulteriori Informazioni



P3 - Gegevensplaatjes, weergaven en bedieningsknoppen
P28 - Inleiding
- Verklaring van illustraties
P29 - Beginnen
P30 - Continu gebruik
P31 - Diagnoseprocedure
P32 - Technische specificaties
P33 - Aanvullende informatie



Σ3 - Χαλκομανίες, ενδείξεις και χειρισμοί
Σ64 - Εισαγωγή
- Οδηγός εικόνων
Σ65 - Πώς αρχίζετε
Σ66 - Συνεχής λειτουργία
Σ67 - Λύση προβλημάτων
Σ68 - Τεχνικές προδιαγραφές
Σ69 - Συμπληρωματικές πληροφορίες



Sida 3 - Dekaler, displayer och reglage
Sida 34 - Inledning
- Bildförklaring
Sida 35 - Förberedelser
Sida 36 - Fortsatt användning
Sida 37 - Felsökning
Sida 38 - Tekniska data
Sida 39 - Tillkommande information

2100 Classic



Introduction

Thank you for choosing the Prestige Medical 2100 *Classic* Autoclave designed to sterilize solid, un-wrapped instruments. After removing from the box, please check for any transit damage. If any damage to the unit is found, please contact your supplier immediately.

Together with this unit and operating manual, you will find the following:

- Electrical mains cord
- Instrument furniture
- Spare fuses x2 (when fitted)
- Performance test certificate
- Certificate of Compliance (230v models only)
- Warranty registration card

Key to Pictures, Diagrams and Symbols

The following descriptions refer to the pictures of the controls, display lights and operating symbols on Page 3 of this manual.

Controls:

- A - Depressurization Valve
- B - Start Cycle Button
- C - Pressure Rise Indicator
- D - Display Panel
- E - Pressure Gauge ("Plus" models only)
- F - Temperature Gauge ("Plus" models only)
- G - WARNING!
The top moulding is not a handle - do not use to remove the lid or lift the autoclave - use the side handles.
HOT PARTS! - Do Not Touch the top moulding during and after a cycle.

Figures:

- F1 - Lid aligned with body, turn clockwise to close
- F2 - Lid in closed position

Warning Symbols:

- W1 - Caution, electric shock hazard
- W2 - Warning, read manual before using the autoclave
- W3 - Warning, unit must be earthed/grounded
- W4 - Warning, heat hazard

Display Lights:

- L1 - Power "On" light - illuminates GREEN
- L2 - Heating Light - illuminates ORANGE
- L3 - Sterilizing light - illuminates YELLOW
- L4 - Sterilizing complete light - illuminates GREEN
- L5 - Fault light - illuminates RED



WARNING!

Do not touch the body or lid as these parts become hot when the autoclave is in operation



Getting Started

Before using the autoclave for the first time please take time to read the following pages to familiarise yourself with the operation of the unit. We strongly recommend that all users of the autoclave are trained in its operation.

The autoclave is very easy to use. By following this simple operating sequence in conjunction with the pictures of the autoclave, it's controls, display panel and operating symbols (page 3), you will be able to ensure your instruments are correctly sterilized every time.

WARNING!

Mains outlet MUST BE EARTHED (GROUNDED)

The mains plug should always be easily accessible as it is to be relied upon as "the means of disconnection"

1. Water

Fill the unit to the water level line on the inside of the chamber with 0.75l of distilled or de-ionised water.

DO NOT USE TAP WATER OR OVERFILL.

2. Loading



Place UNWRAPPED and WASHED instruments ONLY plus a TST Strip into the instrument basket or cassettes, and place in the unit.

If using a basket, place the metal 'V' support in the bottom of the unit so ensuring instruments and TST strips are not in the water. TST strips should be placed as near to the centre of the load as possible.

WARNING!

The manufacturers of instruments should be consulted about their suitability for autoclaving and the maximum temperature which the instruments can withstand.

3. Closing.



Always place the lid on the autoclave with the Depressurization Valve (A) open.

Align the lid and body (F1) and turn in a clockwise direction ensuring the lid is completely closed (F2).

Close the Depressurization Valve (A) so it is aligned with the "O" on the lid.

Never leave in position as shown in F1.

4. Power Connection.

Attach the cable supplied to the rear of the unit and plug into an EARTHED mains electrical socket of the CORRECT voltage.

Lights: L1 illuminates GREEN

5. Starting.

Start the sterilizing cycle by pressing button (B).

Lights: L1 illuminates GREEN
L2 illuminates ORANGE

-as the temperature rises, air will be displaced by steam through the Air Bleed Device located in the lid, until it closes with an audible "click", sealing the unit.

The Pressure Indicator (C) will rise indicating the unit is now pressurised.

-Sterilizing Temperature is reached when:

Lights: L1 illuminates GREEN
L2 flashes ORANGE
L3 illuminates YELLOW

-Sterilizing Cycle is completed when:

The Buzzer sounds.

Lights: L1 illuminates GREEN

L4 illuminates GREEN

Note: L4 remains GREEN until a new cycle is started or the unit is disconnected from the mains power.

6. Depressurizing.

Once the Sterilizing Cycle has been completed, the unit needs to be depressurized and allowed to cool down before the lid and sterilized instruments can be removed.

The time taken to reach the point at which this is safe to do so can be shortened by Manually Depressurizing the unit.

Open the Depressurization Valve (A) by turning slowly in an anti clockwise direction*.

The Pressure Indicator (C) will drop once the steam has been released.

* Warning: There will be a visual and audible release of steam from the rear of the top moulding..

7. Unlocking.

Once the pressure has been released the lid may be unlocked. Ensure the Depressurization Valve (A) is open. Remove the lid by turning in an anti-clockwise direction (F1).

8. Unloading.

Lift off the lid, gently place upside down on a solid work surface and leave to cool. Ensure the Depressurization Valve (A) is in the closed position to avoid damaging it.

The unit has completed a successful cycle if the "spot" on the TST Indicator Strip has completely changed colour from yellow to purple.*

The basket or cassettes containing the sterilized instruments can now be lifted out of the unit.

To avoid damage, replace the lid as described in step "3"

* Please Note: If the "spot" has not completely changed colour, replace with a new TST strip and start a new cycle. If the "spot" fails to change colour for a second time, do not use the unit until it has been checked by a qualified engineer.

DO NOT USE THE INSTRUMENTS IF A COMPLETE STERILIZATION CYCLE HAS NOT BEEN ACHIEVED.

Continued Operation

To ensure your autoclave gives you the years of service for which it was designed, it is important to remember a few "do's" and "don't's" with regards to the operation of the unit and to carry out the simple care and maintenance procedures on a weekly basis.

Do ensure that....

- 1....you read these instructions and always follow the operating sequence.
- 2....the work surface on which you will place the autoclave is flat and solid.
- 3....the instruments are designed to withstand the sterilizing temperatures selected, are thoroughly cleaned and rinsed before sterilizing, and are not any longer than the length, or exceed the load weight, specified - see "Technical data" section.
- 4....the water level is maintained regularly with clean distilled or de-ionised water only.
- 5....the unit is in a "draught free" environment and is positioned not less than 250mm from adjacent walls.
- 6....you only use green sealing gasket (219500) and that it is changed at the end of its life, if visibly damaged, or when  shrinkage has occurred, see "Fault mode - 5".
- 7....the lid is securely closed when the unit is not in use, to avoid the risk of accidental damage.
Never leave in position as shown in F1.
- 8....you quote your model details, serial number and date of purchase when contacting Prestige Medical or your supplier.

Do not....

- 1....touch the unit whilst in operation - it gets HOT.
- 2....attempt to remove the lid during operation.
- 3....lose this operating instruction manual
- 4....add any chemicals whatsoever to the water.
- 5....attempt to sterilize volatile substances toxic materials or inappropriate loads.
- 6....place the unit on heat sensitive surfaces polished wood or glass.
- 7....open the Depressurization Valve (A) during the sterilization cycle.
- 8... leave the Depressurization Valve (A) in the "open" position when placing the lid upside down on a work surface.
- 9... immerse the unit or electrical cord in water when cleaning.
10. use abrasive materials or lubricants when cleaning.
- 11..drop or abuse the unit.
- 12..use in areas of risk associated with flammable materials or gasses.
- 13..attempt to change fuses until the unit has been unplugged from the mains. Only qualified persons should change fuses.
- 14..reach over the unit when removing cover, to do so may cause burns from rising heat and steam.
- 15..press start button once cycle has been started as this will re-set the cycle timer to zero.

Care and Maintenance.

WARNING!
Disconnect the autoclave from the mains power supply before cleaning.

Green Sealing Gasket.

1. Remove from inside the lid and clean with warm, soapy water.
2. Rinse thoroughly, shake dry, do not wipe.
3. Replace in the lid by tucking evenly under all lugs starting at the Gasket Offset Device. It may appear slightly wrinkled until used.
4. Replace gasket when it begins to show signs of leakage.

Autoclave.

6. If a new gasket leaks, or a persistent leak develops, gently clean the sealing surface of both the lid and body of the unit with a "Scotchbrite" scrubbing pad making sure you do not remove any metal. Rinse both surfaces but do not dry.
7. Clean both interior and exterior with warm, soapy water ensuring the electrical parts are kept dry.
8. Monitor the first cycle of the day to check the Air Bleed Device, which is located inside the lid, audibly "clicks" shut.
9. Prestige Medical recommend that your unit is calibrated at six monthly intervals.
10. Lubricate underside of body lugs with "vaseline" if the lid becomes stiff.

DO NOT LUBRICATE GASKET

Troubleshooting

In the event of a fault occurring during any stage of the units operation, identify the fault by referring to the descriptions below. The fault can be rectified by following the Fault Remedy applicable to the problem incurred.

Fault Indication/Description/Remedy

Fault 1: No Power to Unit

Light: L1 fails to illuminate.
Blown fuse / Defective Socket /
Mains not connected

Ensure mains lead is connected.
Check / replace fuse:
Check power to socket:

Fault 2: Low water or Boil dry.

Light: L5 flashes RED

Allow unit to cool before refilling to the correct level.
Disconnect from mains then reconnect and repeat cycle.
If the fault repeats with sufficient water, arrange for a
service engineer to visit.

Fault 3: Sterilization failed to be achieved.

Light: L4 fails to illuminate GREEN and there is no audible
buzzer

Disconnect from mains then reconnect and repeat cycle.
If the fault repeats arrange for a service engineer to visit.

Fault 4: Incomplete Sterilization cycle

TST strip fails to change / completely change colour.

Check expiry date of TST strips.
Disconnect from mains then reconnect and repeat cycle.
If the fault repeats arrange for a service engineer to visit.

Fault 5: Steam or water leaks from under the lid

i) Worn or dirty gasket.

Wash gasket and sealing surfaces on the body and lid as
described under "Care and Maintenance".
If the fault persists, replace with a new gasket.

ii) Incorrectly closed lid.

Ensure the unit is fully depressurized by opening
Depressurization Valve (A). Remove lid and re-fit carefully.
Disconnect from mains, reconnect and repeat cycle.

Fault 6: Excessive steam or water leaking from Depressurization Valve (A).

Depressurization Valve (A) in "OPEN" position.
Close Depressurization Valve (A).

Fault 7: Incorrectly reading gauges.*

* Only applicable to "Plus" models.
Disconnect from mains then reconnect and repeat cycle.
If the fault repeats arrange for a service engineer to visit.

Technical Specifications

In the unlikely event that something should go wrong, we have incorporated a number of safety features to ensure that your autoclave remains safe at all times.

Safety features.

1. Located to the rear of the lid, beneath the cover, is a spring called the Gasket Offset Device (GOD Spring), designed to prevent pressure building up if the lid has been incorrectly fitted.
DO NOT TAMPER WITH THIS SAFETY DEVICE
2. If for any reason, the temperature falls below the minimum required sterilizing temperature, resulting in the Sterilizing Temperature light (L3) switching off, the cycle timer will re-start from zero once the correct temperature has been restored.
3. If there is an electrical or electronic failure resulting in a build up of pressure - in excess of normal operating pressure - one or all of the following safety features will be activated.
 - i) Depressurization Valve (A) will loudly and rapidly "vent" steam.
 - ii) The gasket will "extrude" through the slot in the rear of the lid rapidly releasing excess pressure and steam. 
 - iii) A non - resettable thermal fuse located in the base of the unit will "melt" at a pre-determined temperature, disconnecting the power.

Should any of the devices listed activate, please observe the following steps:

- a) Do not touch the unit.
- b) Switch off at the wall socket and un-plug
- c) Allow temperature and pressure to drop before
 - i) touching the unit
 - ii) removing your instruments
- d) Do not attempt to re-start the unit
- e) Arrange for an immediate service.

Technical Specifications - Standard Body Autoclaves

Height.....	335 mm	Capacity.....	9 litres	Max' Instrument	
Width.....	340 mm	Internal chamber		length.....	228 mm
Net Weight.....	4.5 kilos	Dimensions (d/h).....	210/230 mm	Max' Load Weight.....	3.0 kilos

- Extended Body Autoclaves

Height.....	420 mm	Capacity.....	12 litres	Max' Instrument	
Width.....	340 mm	Internal chamber		length.....	290 mm
Net Weight.....	4.5 kilos	Dimensions (d/h).....	210/270 mm	Max' Load Weight.....	4.0 kilos

- Super Extended Body Autoclaves

Height.....	582 mm	Capacity.....	16.8 litres	Max' Instrument	
Width.....	340 mm	Internal chamber		length.....	460 mm
Net Weight.....	4.5 kilos	Dimensions (d/h).....	210/444 mm	Max' Load Weight.....	4.0 kilos

Fuses - Located under the control module, fuses F10A, 32 x 6.3mm, ceramic sand filled, Mains plug top fuse (User replaceable), F13A to BS1362 UK ONLY.

Rating - Models are rated continuously for intermittent use.

Body - Deep drawn aluminium.

Lid - Aluminium.

Heater - Externally surface mounted mechanically fixed electric element.

Temperature Cut Out - Thermal fuse.

Pressure - Calibrated pressure release valve.

Max. Single Fault Temperature - 133.3°C

Over Voltage Category - Group II

Pollution Degree - Group 2

Environment Conditions - indoor use - temperature 5°C to 40°C - altitude up to 2000m - maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C. - mains supply voltage fluctuations not to exceed + 10% of the nominal voltage

Input Connections - Mains inlet socket 'hot' format conforming to IEC 302.

Safety Shut Down - See 'Temperature Cut Out'.

Packaging - All packaging materials are recyclable.

Additional Information

Spares

Only those spare parts supplied or specified by Prestige Medical should be used in the maintenance of the autoclave. Use of unauthorised parts will invalidate any warranty given and may adversely affect the performance and safety of the unit.

Accessories

A range of accessories are available for your autoclave as described below and pictured on page 76. Contact your supplier for full details.

- 1 - 219294 - Lifting device
- 2 - 219293 - General Instrument Tray
- 3 - 219292 - Standard Basket
- 4 - 219295 - 'V' Support
- 5 - 229310 - Dental Instrument Tray
- 6 - 219291 - Cassette box
- 7 - 219290 - Cassette Rack
- 8 - 219296 - Extended Basket
- 9 - 219500 - Green Silicone Sealing Gasket
- 10 - 219277 - TST Indicator strips
- 11 - 219258 - Cord Set UK
- 12 - 219297 - Cord Set EURO
- 13 - 219299 - Cord Set UL
- 14 - 219509 - Extra large basket (210030/34only)

Warranty

Prestige Medical shall, in the first 12 months from the date of purchase, repair or replace free of charge any parts which prove to be defective in workmanship and/or materials. The heating element (only) is covered by a lifetime guarantee.

Prestige Medical shall not be so liable in the event that the purchaser has failed to adhere to the instructions contained herein or if the autoclave has been abused, interfered with, altered, repaired or serviced by any unauthorised party this may also result in the protection provided by the equipment being impaired.

This warranty excludes the gasket, all internal furniture and consumables.

Consumer's statutory rights are not affected.

Product decontamination.

Should the unit require repair, it must be decontaminated in accordance with a recognised procedure prior to return or on-site repair. A statement of equipment contamination status must be available with the product. (Details of a suitable procedure are available on request).

Approvals.

Note: Please ask for details of approvals specific to your model.

BS3970 Part 4 (pressure vessel only)

CE Mark - Medical Device Directive (93/42/EEC)







Cleaning materials:

- Mild washing up liquid.
- Non-abrasive cream cleaner.
- Disinfectant diluted in water - non chlorine base

Packaging.

All packaging materials used are recyclable, please dispose of accordingly.

PLEASE NOTE: English is the original language for the purposes of these instructions. All other languages are translations from the English text.



- GB** Quote your model details, serial number and date of purchase when contacting Prestige Medical or your supplier
- F** Donnez les informations concernant le modèle, le numéro de série et la date d'achat lorsque vous contactez Prestige Medical ou votre fournisseur local.
- D** Die Modellnummereinheiten, Seriennummer und Kaufdatum bei Kontaktaufnahme mit Prestige Medical oder Ihrer Vertretung bereithalten.
- DK** De angiver modeldetaljer, serienummer og købsdato når De kontakter Prestige Medical eller Deres leverandør.
- NL** De details van uw model, het serienummer en datum van aankoop vermeldt wanneer u contact opneemt met Prestige Medical of uw leverancier.
- S** Ange modellens detaljer, serienummer och inköpsdatum vid kontakt med Prestige Medical eller leverantören.
- N** Oppgir modelldetaljene, serienummeret og kjøpedatoen når du kontakter Prestige Medical eller leverandøren din.
- SF** Mainitset käytössäsi olevan autoklaavin mallia koskevat tiedot, sarjanumeron ja ostopäivän, kun otat yhteyttä -yhtiöön tai laitteen myyjään.
- E** Menciona los detalles del modelo, número de serie y fecha de compra cuando contacte con Prestige Medical o con su proveedor.
- P** Cita os detalhes do modelo, número de série e a data da compra quando contactar Prestige Medical ou o seu fornecedor.
- I** Di indicare i dettagli relativi al modello, il numero seriale e la data di acquisto quando ci si rivolge alla Prestige Medical o al proprio fornitore.
- GR** αναφέρετε τις λεπτομέρειες του μοντέλου σας, τον αύξοντα αριθμό και την ημερομηνία αγοράς όταν επικοινωνείτε με την Prestige Medical ή με τον προμηθευτή σας.

Model no.

Serial no.

Prestige Medical subsidiary / dealer address

